

Iraq Rebuilding Office – Scandinavia i/s

مكتب إعادة أعمار العراق – أسكندينافيا

Trade consultant, sale & troubleshooting

استشارات هندسية / طبية ووكالات تجارية عالمية



E-mail: info@irak-genopbygning.dk

Homesite: www.rebuildingiraq.cjb.net - www.irak-genopbygning.dk

Phone: (+45) 30 22 23 75

Others

This is a selection from our wide range of used types of medical equipments.
The equipments are available from our directly stock and ready to be shipped immediately.

Please do not hesitate to ask for full details and prices.

Welch Allyn Video Colposcope



The Welch Allyn Video Colposcope is the world's first "true" video colposcope. For true video imaging, the Welch Allyn Video Colposcope is a fully self-contained, stand-alone unit that incorporates light, power, and video into a unique, compact system. This breakthrough technology makes all the benefits of high-resolution video imaging more convenient, affordable, and satisfying for both the patient and the physician. Unique to the video system, the Polarization filter removes glare from image for accurate tissue assessment. The patented green filter electronically removes red from image with no light loss. Magnification index displayed on screen. Solarc lamp generates light 50% brighter and whiter than any halogen light source. Unique CCD/optics design 1/4" CCD offers full screen, high-resolution imaging. Compact design with ambidextrous controls. Creates rilliant, high-resolution images because it has no couplers or beam splitters to degrade the image quality or field of view. Solarc illumination provides white, bright light for true color rendition with no light loss at any magnification. Sleek, compact design with ambidextrous controls improves access to the patient. Patent-pending electronic green filter removes red from image with no loss of light. Images can be captured, stored, and printed directly from the instrument. Motorized, digitally controlled camera module provides step-free continuous zoom and focus. Adjustable illumination beam director rotates to eliminate shadows

during any procedure. Simultaneous viewing of video images facilitates group participation among doctors, patients, and staff.

Polarization filter.

Sony 14-inch, high-resolution monitor.

Conforms to UL 2601-1; CAN/CSA C22.2 No. 601.1.

Easy to use.

Ulrich XD5500 CT Contrast Injector



ProAct Workstation

*Software: RadWorks 5.X Version 1.3
3D*

Looking forward to hear from you

Best regards

Iraq Rebuilding Office – Scandinavia